

Manufacturer's Recommendations for Alternate Dental CBCT QA Program

Acteon Group: Model X-Mind Trium

Table 6

Medical Physicist's Computed Tomography QC Survey

Item	Required Test or Procedure	Frequency	Substitute Test or Procedure	Standard
1	Scan Increment Accuracy	Initial & Annually	None – Not Applicable	None – Not Applicable
2	Scan Localization Light Accuracy	Initial & Annually	X-Mind Trium Service Manual 13.1 Adjustment of the Lasers Pages: 139, 144, 145, 150 to 155 (Appendix C) X-Mind Trium User Manual 11.10 Laser Page: 60 (Appendix B)	< 1mm
3	Patient Dose (Multiple Scan Average Dose) MSAD or Computed Tomography Dose Index-CTDI	Initial & Annually	X-Mind Trium User Manual 11.16.2. Dosimetric Indications Pages: 63 to 81 (Appendix B)	See reference dose tables.
4	Pre-Patient Collimation Accuracy	Initial & Annually	X-Mind Trium Service Manual 12.7. Motorized Collimation Calibration Pages: 106 to 113 (Appendix C)	Repeat steps from h to l until the X-ray beam covers correctly the image detector.
5	Contrast Scale	Initial & Annually	X-Mind Trium Quality Assessment Manual 5. CBCT Quality Check 5.2.3. Contrast to Noise Ratio (CNR) 5.2.5. Modular Transfer Function Pages: 5 to 21, 23 to 25, 27 (Appendix A)	Trium: Man-Medium, Medium Quality, FOV 80x80 OR Trium TLD: Man-Medium, Medium Quality, FOV 80x90 CNR > 5.0 $V_{10\%} [LP/mm] > 1.6$ $V_{50\%} [LP/mm] > 0.5$ Trium: Man-Medium, Medium Quality, FOV 110x80 OR Trium TLD: Man-Medium, Medium Quality, FOV 110x90

				<p>CNR > 4.0</p> <p>$V_{10\%}$ [LP/mm]: > 1.5</p> <p>$V_{50\%}$ [LP/mm]: > 0.4</p>
6	CT Number for Water	Initial & Annually	<p>X-Mind Trium Quality Assessment Manual</p> <p>5. CBCT Quality Check</p> <p>5.2.4 Homogeneity</p> <p>Pages: 5 to 19, 22, 23 & 27 (Appendix A)</p>	<p>The software computes automatically the homogeneity value.</p> <p>Trium: Man-Medium, Medium Quality, FOV 80x80 OR Trium TLD: Man-Medium, Medium Quality, FOV 80x90</p> <p>Homogeneity > 12</p> <p>Trium: Man-Medium, Medium Quality, FOV 110x80 OR Trium TLD: Man-Medium, Medium Quality, FOV 110x90</p> <p>Homogeneity > 10</p>
7	Slice Thickness	Initial & Annually	None – Not Applicable	None – Not Applicable
8	Field Uniformity	Initial & Annually	<p>X-Mind Trium Quality Assessment Manual</p> <p>5. CBCT Quality Check</p> <p>5.2.4 Homogeneity</p> <p>Pages: 5 to 19, 22, 23 & 27 (Appendix A)</p>	<p>The software computes automatically the homogeneity value.</p> <p>Trium: Man-Medium, Medium Quality, FOV 80x80 OR Trium TLD: Man-Medium, Medium Quality, FOV 80x90</p> <p>Homogeneity > 12</p> <p>Trium: Man-Medium, Medium Quality, FOV 110x80 OR Trium TLD: Man-Medium, Medium Quality, FOV 110x90</p> <p>Homogeneity > 10</p>
9	Low Contrast Resolution	Initial & Annually	<p>X-Mind Trium Quality Assessment Manual</p> <p>5. CBCT Quality Check</p> <p>5.2.3. Contrast to Noise Ratio (CNR)</p> <p>Pages: 5 to 21, 27 (Appendix A)</p>	<p>Trium: Man-Medium, Medium Quality, FOV 80x80 OR Trium TLD: Man-Medium, Medium Quality, FOV 80x90</p> <p>CNR > 5.0</p>

				Trium: Man-Medium, Medium Quality, FOV 110x80 OR Trium TLD: Man-Medium, Medium Quality, FOV 110x90 CNR > 4.0
10	High Contrast Resolution	Initial & Annually	X-Mind Trium Quality Assessment Manual 5. CBCT Quality Check 5.2.5. Modular Transfer Function Pages: 5 to 19, 23 to 25, 27 (Appendix A)	Trium: Man-Medium, Medium Quality, FOV 80x80 OR Trium TLD: Man-Medium, Medium Quality, FOV 80x90 $V_{10\%} \text{ [LP/mm]} > 1.6$ $V_{50\%} \text{ [LP/mm]} > 0.5$ Trium: Man-Medium, Medium Quality, FOV 110x80 OR Trium TLD: Man-Medium, Medium Quality, FOV 110x90 $V_{10\%} \text{ [LP/mm]} > 1.5$ $V_{50\%} \text{ [LP/mm]} > 0.4$
11	Noise	Initial & Annually	X-Mind Trium Quality Assessment Manual 5. CBCT Quality Check 5.2.3. Contrast to Noise Ratio (CNR) Pages: 5 to 21, 27 (Appendix A)	Trium: Man-Medium, Medium Quality, FOV 80x80 OR Trium TLD: Man-Medium, Medium Quality, FOV 80x90 CNR > 5.0 Trium: Man-Medium, Medium Quality, FOV 110x80 OR Trium TLD: Man-Medium, Medium Quality, FOV 110x90 CNR > 4.0
12	Scan Protocol Review	Initial & Annually	See N.J.A.C. 7:28-22.10(a)	N.J.A.C. 7:28-22.10(a)
13	Review of Facility and Technologist's QC Tests	Initial & Annually	See N.J.A.C. 7:28-22.10(a)	N.J.A.C. 7:28-22.10(a)
14	Physicist Report and Recommendations	Initial & Annually	See N.J.A.C. 7:28-22.10(a)	N.J.A.C. 7:28-22.10(a)